

D-KAT Digital Keeler Applanation Tonometer



D-KAT

The Digital Keeler Applanation Tonometer (D-KAT) is intended for measuring intraocular pressure to aid in the screening and diagnosis of glaucoma.

This device is intended to be used only by suitably trained and authorised healthcare professionals.



Key features

- British built
- LED display for ease of use in darkened environments
- Lightweight and portable
- Display reading to 1 decimal point
- Low battery usage
- Less moving parts, meaning increased reliability and less maintenance

Digital Keeler Applanation Tonometer



Specifications

D-KAT Installation

| KEELER DIGITAL T-TYPE: for H-Style illumination systems | Fitted to the guide plate on the optical axis for the microscope and illumination unit arm |
|--|--|
| KEELER DIGITAL R-TYPE: for some Z-Style and H-Style illumination systems | Mounted on tonometer post or microscope |
| Measurement | |
| MEASUREMENT RANGE | 5 - 65 mmHg |
| STANDARD DIVERGENCE | 0.49 mN or 1.5% of measurement value, whichever is the greater |
| Operation | |
| OPERATING TEMPERATURE RANGE | From 10°C to 35°C |
| MEASUREMENT ACCURACY | ≤ 0.49 mN |
| Net Weight | |
| KEELER DIGITAL T-TYPE | 0.37 kg (without accessories) |
| KEELER DIGITAL R-TYPE | 0.58 kg (without accessories) |
| | |
| D-КАТ Т-Туре | |



Clear LED display

Part Numbers

Keeler D-KAT Tonometers

| 2414-P-2032 | Keeler D-KAT T-Type |
|-------------|---------------------|
| 2414-P-2042 | Keeler D-KAT R-Type |

Accessories

| 2414-P-8037 | Tonometer doubling prism |
|-------------|--------------------------|
| 2414-P-5005 | Calibration arm assembly |
| 2414-P-5032 | T-Type guide plate |
| 2414-P-5042 | R-Type post |
| 1909-P-7129 | AA Battery |
| | |





Doubling Prism

Calibration arm assembly



T-Type guide plate

R-Type post

D-KAT Overview

D-KAT

The Digital Keeler Applanation Tonometer (D-KAT) is a screening device used to measure intraocular pressure which is one of the key biomarkers considered in diagnosing glaucoma.

D-KA1

The D-KAT is a Digital Goldmann Tonometer which uses the same principle found in non-digital tonometers, however removes the potential error reading from the dial indicating the IOP.

The D-KAT comes in two variants, the R-Type and the T-Type, both of which can be used on tower style Slit Lamps.

What's in the box?

- D-KAT
- Instructions for use USB
- AA battery
- Calibration key
- Glide plate (T-Type only)
- Tonometer post (R-Type only)
- Tonometer mount screw (R-Type only)
- Reusable tonometer prism

More Keeler Instruments

With over 100 years of industry experience, Keeler has become a household name among optometrists and ophthalmologists worldwide.

View our full range at www.keeler.co.uk | www.keelerusa.com

02012021 REV A

Got a question? - Call us

UK: +44 (0)1753 857177 India: +91 99 3031 1090

USA: +1 610 353 4350 Brazil: +55 11 4302 6053

China: +86-18512119109





Manufacturing address 🖬 Keeler Ltd. Clewer Hill Road Berkshire SL4 4AA UK

CE₁₆₃₉