

**Z-KAT**Digital Keeler Applanation Tonometer



# Z-KAT

The Digital Keeler Applanation Tonometer Z-type (Z-KAT) is intended to determine the intraocular pressure (IOP) by exerting an external force against the eye which provides reading of the resistance of tunica of the eye of deformation.



### Key features

- British built
- LED display for ease of use in darkened environments
- Lightweight and portable
- Display reading to 1 decimal point
- Low battery usage
- Less moving parts, meaning increased reliability and less maintenance

# Digital Keeler Applanation Tonometer



## **Specifications**

### **Z-KAT Installation**

KEELER DIGITAL Z-TYPE for Z-Style illumination systems

Mountable on post on microscope

INSTALLATION

Assembled on the tonometer post on the slit lamp

#### Measurement

MEASUREMENT RANGE

5 - 65 mmHg

STANDARD DIVERGENCE

0.49 mN or 1.5% of measurement value, whichever is the greater

From 10°C to 35°C

#### Operation

OPERATING

TEMPERATURE RANGE

 $\begin{array}{ll} \textbf{MEASUREMENT ACCURACY} & \leq 0.49 \ mN \end{array}$ 

**Net Weight** 

KEELER DIGITAL Z-TYPE

0.6 kg (without accessories)





Clear LED display

### Part Numbers

#### **Keeler Z-KAT Tonometer**

2414-P-2010 Keeler Z-KAT Tonometer

#### Accessories

2414-P-8037 Tonometer doubling prism
2414-P-5045 Calibration arm assembly
2414-P-5073 Z-type post
1909-P-7129 AA Battery







Calibration arm assembly



Z-Type post



### **Z-KAT Overview**

The Keeler Digital Applanation Tonometer Z-style (Z-KAT) is a screening device used to measure intraocular pressure which is one of the key biomarkers considered in diagnosing glaucoma.

The Z-KAT is a digital Goldmann Tonometer which uses the same principal found in non-digital tonometers, however removes the potential error reading from the dial indicating the IOP.

# What's in the box?

- Z-KAT
- Instructions for use
- AA battery
- Calibration key
- Tonometer post
- Tonometer mount screw
- Reusable tonometer prism
- Digital spirit level
- Allen key

### More Keeler Instruments

With over 100 years of industry experience, Keeler has become a household name among optometrists and ophthalmologists worldwide.

View our full range at www.keeler.co.uk | www.keelerusa.com

### Got a question? - Call us

UK: +44 (0)1753 857177 India: +91 99 3031 1090

USA: +1 610 353 4350 Brazil: +55 11 4302 6053

China: +86-18512119109





**C**€ 1639



